

DISCOUNTED PRICING OFFERS SAVINGS AND VALUE

CLASSIC ELECTROPHYSIOLOGY DESIGNS AVAILABLE THROUGH A SINGLE SYSTEM CONFIURATION

EACH BUNDLED SYSTEM COMES WITH 2 MANIPULATORS

EASY TOGGLE SELECTION OF ACTIVE COMPONENT

SINGLE ROE INPUT DEVICE CONTROLS MANIPULATORS AND MOTORIZED COMPONENTS PROVIDING EASE OF OPERATION

ALL FEATURES OF MANIPULATORS, STAGES AND TRANSLATORS ARE RETAINED

SYSTEMS WITH PLATFORMS INCLUDE ROTATION BASES FOR MOUNTING MANIPULATORS

SIMPLE USB INTERFACE

Shown: 1078-385-N65. Headstages and microscope not included)

ELECTROPHYSIOLOGY BUNDLED SYSTEMS

MANIPULATOR / STAGE / PLATFORM BUNDLED SYSTEMS

Over the last 20 years, electrophysiology rigs have developed spontaneously around several platforms. In all, the central concept is that manipulators and the sample are fixed in one frame and the microscope is fixed in a second frame. These two separate frames of reference are then moved with respect to each other in order to locate cells. Sutter Instrument is now offering several systems or "Big Kits" that include dual manipulators and either manipulator stands and scope translator or large moving stage and dual manipulators systems as one big kit. By specifying and offering the Big Kit systems, Sutter is able to extend significant price savings over the same components purchased separately.

Sutter San Francisco... the classic moving scope design with manual translator and three columns and two manipulators. Developed by labs at UCSF and later used throughout the Bay area, the U.S., and the world. The time-honored method allows for easy hand access to controls on the scope and manipulators for pipette exchange and other adjustments.

Sutter Alcatraz... another classic design. Relies on a manual scope translator but uses a single large fixed plate stage rather than separate stands for manipulators. Many electrophysiologists want the ability to add additional manipulators or other devices to a single large platform. Both the San Francisco and the Alcatraz are available with a motorized (rather than manual) translator.

Sutter Long Island ... intended specifically for multipipette electrophysiology on two-photon microscopes. Designed in conjunction with Northwestern University and first used for the CSHL Imaging course. The large moving stage can easily hold an in-vivo or slice preparation and two or more manipulators. The scope is allowed to be fixed to the table top which is required for scopes attached to two-photon sources. This same design is also good for systems where the scope has become too large to move easily like spinning disc confocals and large camera systems. The stage and manipulator systems have been used on many two-photon and other scopes including those from Intellegent Imaging Innovations, Zeiss, Olympus and Leica. This configuration is also used by researchers in custom built two-photon microscopes.

Software compatibility. Motorized translation and motorized stage systems are compatible with a wide range of commercial and freeware imaging software platforms including Intelligent Imaging Innovations' Slidebook, Scan Image, MicroManager and Sutter's **MCS**, **MOM**[®] computer system. Please contact Sutter for a current listing of software suites that are compatible with our products.

SAN FRANCISCO

MT-1000 with manual X-Y translator, micromanipulators 1000-325-(scope suffix) With two MP-225 manipulators 1000-385-(scope suffix) With two MP-285 manipulators

MT-2200 with motorized X-Y translator, micromanipulators2200-325-(scope suffix)With two MP-225 manipulators2200-385-(scope suffix)With two MP-285 manipulators

ALCATRAZ

MT-1078 with manual X-Y translator, micromanipulators 1078-325-(scope suffix)¹ With two MP-225 manipulators 1078-385-(scope suffix)¹ With two MP-285 manipulators

MT-2278 with motorized X-Y translator, micromanipulators 2278-325-(scope suffix)[†] With two MP-225 manipulators 2278-385-(scope suffix)[†] With two MP-285 manipulators

LONG ISLAND

MT-78 with motorized X-Y	' translator, micromanipulators
78-325 ¹	With two MP-225 manipulators
78-385 ¹	With two MP-285 manipulators

MICROSCOPE MODELS -

Use this	s scope suffix when (ordering	
Z25	Zeiss Axioskop 2 FS	Z25	Leica DMLFS
Z45	Zeiss Axio Examiner	Z25	Leica DM6000 FS
Y51	Olympus BX51WI		DM6 FS
Y51-FD	Olympus BX51WI	Z25	Nikon FN1
	(includes focus drive)		
Y53	Olympus BX53/63		

¹ Specify stage insert type when ordering.



SUTTER INSTRUMENT

One Digital Drive • Novato • CA 94949 • Phone 415.883.0128 Fax 415.883.0572 • Web www.sutter.com • Email info@sutter.com



SUTTER INSTRUMENT

Others may be available. Please call Sutter for details.

Z25 Zeiss Axioskop 2FS
Z45 Zeiss Axio Examiner
Y51 Olympus BX51WI
Y51-FD Olympus BX51WI (includes focus drive)
L30 Leica DM1FS
L35 Leica DM6000FS, Leica DM6 FS
N65 Nikon FNI

MICROSCOPE SUFFIXES