



## SUTTER INSTRUMENT MICROMANIPULATORS

Features	MPC-325	MPC-385	MPC-365	MP-285
Manipulators per Controller	2	2	2	1
Manipulators per ROE	4	4	4	1
Memory Moves	2 (HOME & WORK)	2 (HOME & WORK)	2 (HOME & WORK)	500 & HOME
Travel	25 mm on all three axes	25 mm on all three axes	25 mm on X and Z axes 12.5 mm on Y axis	25 mm on all three axes
Interface	USB	USB	USB	Serial
Expandable	Yes	Yes	Yes	No
Step Size	62.5 nm	62.5 nm	62.5 nm	40 nm and 200 nm
Selectable Speeds	6	6	6	14
Max. Speed	3 mm/sec	5 mm/sec	5 mm/sec	2.9 mm/sec
Type of Bearing	Ball Bearing	Cross Roller	Cross Roller	Cross Roller

**SUTTER INSTRUMENT**<sup>®</sup>

One Digital Drive • Novato • CA 94949 • Phone +1.415.883.0128  
 Fax +1.415.883.0572 • Web [www.sutter.com](http://www.sutter.com) • Email [info@sutter.com](mailto:info@sutter.com)