



SUTTER INSTRUMENT MICROPIPETTE PULLERS

Features	P-97	P-2000 (G/F)	P-30
Heat Source	Platinum Filament	CO ₂ Laser	Platinum or Nichrome Filament
Glass Type	Borosilicate & Aluminosilicate	Quartz Glass or Fiber* Borosilicate & Aluminosilicate	Borosilicate & Aluminosilicate
Max. Glass Size (OD)	4mm	1.5mm Quartz, 1.8mm Borosilicate 0.6mm Fiber/Capillary*	2mm
Tip Size	0.06μ - 3μ	0.03μ - 5μ	0.3μ - 2μ
Max. Taper Length	2cm	1.8cm	1cm
Number of Programs	100	100	NA
Program Lines	8	8	2 Stages (Manual)
Type of Cooling	Compressed Dry Air w/Humidity Control Chamber	NA (Laser On/Off)	None
Usable Pipettes/Pull	2	2	1
PRICE	\$6,700	\$13,500	\$3,250

*Applicable only for the P2000/F

SUTTER INSTRUMENT

One Digital Drive • Novato • CA 94949 • Phone 415.883.0128
Fax 415.883.0572 • Web www.sutter.com • Email info@sutter.com