



## SUTTER INSTRUMENT MICROPIPETTE PULLERS

Features	P-1000	P-97	P-2000 (G/F)	P-30
<b>Heat Source</b>	Platinum Filament with Safe Heat Mode	Platinum Filament	CO <sub>2</sub> Laser	Platinum or Nichrome Filament
<b>Glass Type</b>	Borosilicate and Aluminosilicate	Borosilicate and Aluminosilicate	Quartz Glass or Fiber <sup>1</sup> Borosilicate and Aluminosilicate	Borosilicate and Aluminosilicate
<b>Max. Glass Size (OD)</b>	4mm	4mm	1.5mm Quartz, 1.8mm Borosilicate 0.6mm Fiber/Capillary <sup>1</sup>	2mm
<b>Tip Size</b>	0.06μ - 3μ	0.06μ - 3μ	0.03μ - 5μ	0.3μ - 2μ
<b>Max. Taper Length</b>	2cm	2cm	1.8cm	1cm
<b>Number of Programs</b>	100	100	100	NA
<b>Program Lines</b>	4 plus Line Repeat Mode	8	8	1 Stage or 2 Stage with manual adjustment
<b>Type of Cooling</b>	Compressed Dry Air with Humidity Control Chamber	Compressed Dry Air with Humidity Control Chamber	NA (Laser On/Off)	None
<b>Usable Pipettes/Pull</b>	2	2	2	1
<b>PRICE<sup>2</sup></b>	\$7,950	\$7,590	\$14,595 (glass) \$15,595 (fiber)	\$3,519

1. Applicable only for the P2000/F

2. US prices. International prices are 5% higher.

**SUTTER INSTRUMENT**

One Digital Drive • Novato • CA 94949 • Phone 415.883.0128  
Fax 415.883.0572 • Web [www.sutter.com](http://www.sutter.com) • Email [info@sutter.com](mailto:info@sutter.com)