



CE DECLARATION OF CONFORMITY

Application of Council Directives:
2004 / 108 / EC (EMC) and 2006 / 95 / EC (LVD)

Manufacturer's Name: Sutter Instrument Company

Manufacturer's Address: One Digital Drive
Novato, CA. 94949 USA
Tel: +1 415 883 0128

Equipment Tested: Lambda 10-B Optical Filter Switching & Shuttering System

Model(s): Controller: LB10-B/IQ
10-Position Filter Wheels: LB10-NW (25mm (no shutter))
LB10-NWS (25mm, with standard shutter)
LB10-NWIQ (25mm, with SmartShutter)
LB10-W32IQ (32mm, with SmartShutter)


Conforms to Standards: EMC: EN 61326-1:2006, including:
EN 55011:2009 Class A; EN 61000-3-2:2006 + A2:2009, &
EN 61000-3-3:2008
EN 61000-4-2:2009, BS EN 61000-4-3:2006 + A2:2010,
BS EN 61000-4-4:2004 + A1:2010, EN 61000-4-5:2006,
EN 61000-4-6:2009, EN 61000-4-8:2010, &
EN 61000-4-11:2004
LVD (Safety): EN 61010-1:2010

Tested By: ITC Engineering Services, Inc.
9959 Calaveras Road, PO Box 543
Sunol, CA 94586-0543 USA
Tel. +1 925 862 2944 Fax: +1 925 862 9013
Email: itcemc@itcemc.com Web: www.itcemc.com

Year Tested: 2012

Sutter Instrument hereby declares that the equipment specified above was tested and conforms to the Directives and Standards listed above.

Project Engineer:


Jack Belgum, Ph.D.
Senior Vice President

SUTTER INSTRUMENT

ONE DIGITAL DRIVE, NOVATO, CA 94949 USA PHONE: +1 415 883 0128 FAX: +1 415 883 0572
EMAIL: INFO@SUTTER.COM WEB: [HTTP://WWW.SUTTER.COM](http://WWW.SUTTER.COM)