



## SUTTER INSTRUMENT MICROMANIPULATORS

Features	MPC-325	MPC-385	MPC-365	MP-285
<b>Manipulators per Controller</b>	2	2	2	1
<b>Manipulators per ROE</b>	4	4	4	1
<b>Memory Moves</b>	2 (Home & Work)	2 (Home & Work)	2 (Home & Work)	500 & Home
<b>Travel</b>	25mm on all three axes	25mm on all three axes	25mm on X and Z axes 12.5mm on Y axis	25mm on all three axes
<b>Interface</b>	USB	USB	USB	Serial
<b>Expandable</b>	Yes	Yes	Yes	No
<b>Step Size</b>	62.5nm	62.5nm	62.5nm	40nm and 200nm
<b>Selectable Speeds</b>	6	6	6	14
<b>Max. Speed</b>	3mm/sec	5mm/sec	5mm/sec	2.9mm/sec
<b>Type of Bearing</b>	Ball Bearing	Cross Roller	Cross Roller	Cross Roller
<b>PRICE<sup>1</sup></b>	\$7,830 (MPC-325) \$12,505 (MPC-325-2)	\$9,375 (MPC-385) \$15,330 (MPC-385-2)	\$10,150 (MPC-365) \$16,105 (MPC-365-2)	\$9,720

1. US prices. International prices vary by country.

**SUTTER INSTRUMENT**<sup>®</sup>

One Digital Drive • Novato • CA 94949 • Phone: 415.883.0128  
Fax: 415.883.0572 • Web: [www.sutter.com](http://www.sutter.com) • Email: [info@sutter.com](mailto:info@sutter.com)