



**50mm Round Chamber**



Used with following Sutter Instrument Stages

<a href="#">MP-78</a>	<a href="#">MPC-78</a>	<a href="#">MT-78-FS</a>
<a href="#">MT-1078</a>	<a href="#">MT-2078</a>	<a href="#">MT-2278</a>

*Click on product above to view detailed information*

Part No	Material	Overall	ID	OD
X040510	Aluminum	101.6x152.4mm(4x6in)	41.0mm(1.6in)	50.0mm(2.0in)
X040512	Stainless St	101.6x152.4mm(4x6in)	41.0mm(1.6in)	50.0mm(2.0in)

